




# CETONI

## **CE** NEMESYS OEM

### Hardware Manual



ORIGINAL INSTRUCTIONS 2.05 – JUNE 2019



CETONI GmbH  
Wiesenring 6  
07554 Korbussen  
Germany

**T** +49 (0) 36602 338-0

**F** +49 (0) 36602 338-11

**E** [info@cetoni.de](mailto:info@cetoni.de)

[www.cetoni.de](http://www.cetoni.de)

The information and data contained in this document are subject to change without notice. CETONI GmbH is constantly striving to develop all its products. This means that there may be changes in form, equipment and technology. Claims can therefore not be made on the basis of information, illustration or descriptions in these instructions. The description for the product specification in this manual does not constitute an integral part of the contract.


If you control the products with a software from CETONI GmbH, you agree to the applicable license agreement, which can be read in the corresponding software manual. This and all other current product manuals can be found at <https://www.cetoni.de/en/downloads/manuals>.

CETONI GmbH grants its customers the right to reproduce this manual for the purpose of providing technical information to potential users of CETONI products. Extracts of these documents may only be reproduced or transmitted with a precise reference to the authorship of CETONI GmbH, irrespective of the manner in which they are reproduced or transmitted, either electronically or mechanically. Any duplication or use of excerpts for other purposes requires the written permission of CETONI GmbH.

We are always open to comments, corrections and requests. Please send them to [info@cetoni.de](mailto:info@cetoni.de).

The general terms and conditions of CETONI GmbH shall apply. Alternative agreements must be in written form.

Copyright © CETONI GmbH – Automation and Microsystems. All rights reserved.



# 1 Overviews and Indexes

## 1.1 Content

1	Overviews and Indexes	5
1.1	Content	5
1.2	Revision History	8
2	Introduction	9
2.1	Preface	9
2.2	Symbols and Keywords	9
3	Basic Information	10
3.1	Standards and Guidelines	10
3.2	Intended Use	10
3.2.1	General Description of the Component	10
3.2.2	Proper Use	10
3.2.3	Reasonably Foreseeable Misuse	10
3.2.4	Safety Measures	11
3.2.5	Measures for Safe Setup	12
3.2.6	Safe Operation Measures	12
3.2.7	Condition of Devices	13
3.3	Warranty and Liability	13
4	Technical Data	14
4.1	Performance Data	14
4.1.1	Mechanical Data	14
4.1.2	Electrical Data	14
4.1.3	Inputs	14
4.1.4	Outputs	15
4.1.5	Interfaces	15
4.1.6	Ambient Conditions	15

4.1.7	Connections	15
4.2	Dimensional Drawing of Pumps	16
4.3	Standard Valve	17
4.4	Dimensional Drawing of Valve	17
4.5	Dosing Performance	18
5	Transport and Storage	19
5.1	Transport	19
5.2	Maintenance and Care	19
6	Hardware Operation	20
6.1	Installing a Syringe	20
6.2	Fluidic / Valve	23
6.3	Removing Valves	24
7	Electrical Interfaces	26
7.1	Overview	26
7.2	Power Supply (J1 / J2)	27
7.2.1	Pin Assignment	27
7.2.2	Technical Data	27
7.3	CAN Interface (J3 / J4)	28
7.3.1	Pin Assignment	28
7.3.2	Technical Data	28
7.3.3	Connecting neMESYS to CAN bus line CiA DS-102	28
7.4	RS-232 Connection (J5)	29
7.4.1	Pin Assignment	29
7.4.2	Technical Data	29
7.4.3	Connecting neMESYS to a PC	29
7.5	Signal Connection (J6)	30
7.5.1	Pin Assignment	30
7.5.2	Analog Input 1 and 2 (Pin 1 and 2)	31
7.5.3	Digital Inputs 1 and 2 (Pin 3 and 4)	32
7.5.4	Digital Input 3 (Pin 5)	32

7.5.5	Digital Output 1 and 2 (Pin 6 and 7)	33
7.5.6	Digital Output 3	35
7.6	CAN-Bus Termination	36
7.7	Switching the Integrated Valve	36
8	Cable Set	38
8.1	Power Cable	38
8.2	Power Connection Cable	38
8.3	CAN Cabel (terminated DSub Socket)	39
8.4	CAN Cabel (DSub Plug)	40
8.5	CAN Connection Cable	41
8.6	CAN Terminator	42
8.7	RS232 Cable (DSub Socket)	43
8.8	USB-to-CAN Adapter	44
9	Setup & Cable Connection	45
9.1	Mechanical Setup	45
9.2	Introduction to Cable Connection	45
9.3	Step 1 – Connecting the Power Cable	45
9.4	Step 2 – Data Connection	45
9.4.1	Introduction	45
9.4.2	Data Connection Through CAN	46
9.4.3	Data Connection Through RS232	46
9.5	Step 3 – Power Supply of Additional Modules	46
9.6	Step 4 – Data Connection to Additional Modules	46
9.7	Step 5 – CAN Bus Termination	47
10	Disposal	48

## 1.2 Revision History

REV	DATE	MODIFICATION
1.03	17.08.2010	Creation of Manual
2.00	09.01.2015	Revision for new Hardware
2.01	30.03.2015	RS232 cable modified (Dsub jack instead of plug) Added logic table for valve switching
2.02	11.03.2016	New corporate design
2.03	15.11.2017	correction of RS232 cable pin assignment
2.04	17.05.2018	Dimensional drawing updated
2.05	03.06.2018	Fixed inverted pictures in section 8.6 Created navigable table of contents for PDF



# 2 Introduction

## 2.1 Preface

Thank you for purchasing a CETONI product. With this manual we would like to support you in using the device. For additional questions or comments please feel free to contact us directly.

## 2.2 Symbols and Keywords

This manual uses the following symbols, intended to help you navigate the document:



**HINT.** User tips and useful information to simplify the use of the software.



**IMPORTANT.** Important information and additional, particularly useful information.

This symbol does not refer to dangerous or hazardous situations.



**ATTENTION.** Indicates a potentially harmful situation. If it is not avoided, the product or something in its environment may be damaged.



**CAUTION.** Indicates a potentially dangerous situation. If it is not avoided, slight or minor injuries and property damage may result.

# 3 Basic Information

## 3.1 Standards and Guidelines

neMESYS OEM-modules conform to the basic health and safety requirements stipulated by the 2006/42/EC machine directive. The modules are incomplete machines according to the 2006/42/EC machine directive and intended for installation into a machine. According to the directive they do not bear the CE mark.

Integrating the modules into a machine may cause additional risks. It is imperative that you conduct a risk evaluation of your machine with the integrated modules and take any safety measures necessary.



**IMPORTANT.** Please note that you may not use neMESYS OEM-modules until the machine or plant into which the modules are integrated complies with the 2006/42/EC machine directive and the declaration of conformity in attachment II A has been issued.

## 3.2 Intended Use

### 3.2.1 General Description of the Component

neMESYS devices are syringe pumps. They enable emptying and filling of syringes through a linear relative movement between a syringe holder and a piston holder.

The OEM-module is not a complete machine as defined by the machine directive, but an incomplete machine intended for installation into another machine or plant. Please refer to the information provided in chapter 9 of this manual with respect to installation and cable connection.

### 3.2.2 Proper Use

The neMESYS syringe pump system is intended for high-precision and pulsation-free dosing of fluids in a range from nanoliters to milliliters per second. Depending on the device being used, pressures in excess of several hundred bar may be reached. The devices are typically used in a lab-like environment.

### 3.2.3 Reasonably Foreseeable Misuse

Using the devices in applications other than the ones intended for them may create dangerous situations and must be avoided.



**CAUTION.** The devices must not be used as medical devices or for medical purposes.

## 3.2.4 Safety Measures

Operator safety and failure-free operation of the devices can only be guaranteed when using original equipment parts. Only original accessories may be used. Warranty claims are void if damage was caused by using third-party equipment or third-party material.

The devices were developed and designed in such way as to largely rule out dangers, if used properly. Nevertheless, you should observe the following safety measures to rule out any residual danger.

- CETONI GmbH would like to point out the operator's responsibilities when using the device. Local laws and regulations relevant to using this device must be observed. In the interest of a safe work process the operator and user of the device is responsible for observing all relevant laws and regulations.
- The devices must not be used as medical products or for medical purposes.
- The device is designed and approved for operation in systems governed by article 4 section 3 or the pressurized devices directive 2014/68/EU. This means that the devices may not exceed a maximum volume of 1 liter. When using group 1 fluids according to article 13 of the pressurized systems directive 2014/68/EU, the maximum allowable system pressure is 200 bar. It is 1000 bar for group 2 fluids. If the "Technical Data" section stipulates different product-specific maximum pressure values, such values must be observed. Particular attention must be paid to not exceeding the pressure range of the respective sensor. The information provided in the "Technical Data" section with respect to the maximum operating temperature must be observed. CETONI GmbH is not liable for potential impacts caused by the user extending the system with peripheral devices in a way that causes one or both of these values to be exceeded. The operator is responsible for becoming acquainted with the aforementioned pressurized systems directive and observing applicable stipulations.
- Before each operation of the device the user has to ensure that the device is functioning safely and is in proper condition.
- The user must be familiar with the operation of the device and the software.
- Before starting operation, the devices and wiring must be checked for damage. Damaged wires and plugs must be replaced immediately.
- Cables must be routed in such way as to rule out any trip hazards.

- Do not touch moving parts on the devices during operation. There is a risk of crushing.
- Operating the devices in an explosive environment or with explosive substances is prohibited!
- Make sure to wear safety goggles during installation work on the device or when you work with corrosive, hot or otherwise dangerous substances.
- Transport, storage or operation at temperatures below 0°C with water in the fluid channels may cause damage to the devices.

## 3.2.5 Measures for Safe Setup

- Mechanical and electrical installation of pumps may only be done by skilled personnel.
- Make sure that all devices are installed in line with local laws and regulations.
- Please note that, principally, the OEM-pumps cannot be considered fail-safe. You have to ensure that a failing pump will not cause damage to your application by employing suitable safety and monitoring measures.
- Please note that you are not authorized to repair CETONI devices or components.
- Make sure that the power supply is not active and cannot be activated while conducting installation or wiring work on the devices.

## 3.2.6 Safe Operation Measures

### 3.2.6.1 ELECTRONIC EMISSIONS

The neMESYS syringe pump system is designed for usage in any facility directly connected to a public supply network that also supplies buildings used for residential purposes.

### 3.2.6.2 ESD- DISCHARGE

Floors should be wood or concrete or covered with ceramic tiles. If floors are covered with synthetic material, relative air humidity must be at least 30%.

### 3.2.6.3 ELECTRICAL DISTURBANCES

Supply voltage quality should be equal to a typical business or hospital environment.

### 3.2.6.4 MAGNETIC DISTURBANCES

Power lines, including those of other devices, should not be placed near the device or its cables. Mobile two-way radios should be kept at the minimum safe distance from the device and its wiring.

### 3.2.7 Condition of Devices

Despite flawless workmanship, the devices may be damaged during operation. Therefore, you should visually inspect the device components mentioned before each use. Pay particular attention to crushed cables, damaged hoses and deformed plugs. If you find any damage, please refrain from using the devices and contact CETONI GmbH immediately. We will repair your devices as soon as possible. Never attempt to repair a device yourself.

## 3.3 Warranty and Liability

The devices left our facilities in perfect condition and may only be opened by CETONI GmbH. If a device is opened by an unauthorized person, all warranty and liability claims shall be void, in particular those referring to personal injury.

The warranty period is 1 year from the day of delivery. Any work done on the devices within this period shall not extend or renew the warranty.

CETONI GmbH assumes responsibility for its devices with respect to safety, reliability and function only if installation, readjustment, changes, extensions and repairs are done by CETONI GmbH or an authorized party, and if the devices are used in accordance with the user manual.

The neMESYS syringe pump system complies with the applicable safety rules and standards. CETONI GmbH reserves all property rights for the relevant wiring, processes, names, software and devices.

# 4 Technical Data

## 4.1 Performance Data

### 4.1.1 Mechanical Data

<b>DIMENSIONS (L X B X H)</b>	310 x 47 x 130 mm
<b>WEIGHT</b>	≈1100 g
<b>ATTACHMENTS</b>	M3 screws

### 4.1.2 Electrical Data

<b>SUPPLY VOLTAGE VCC (RIPPLE &lt; 10%)</b>	24VDC
<b>PEAK CURRENT DRAW AT 24 VDC</b>	0.6 A
<b>TYPICAL CURRENT DRAW AT 24 VDC</b>	0.3 A

### 4.1.3 Inputs

Analog input 1 (AI1)	resolution 11-bit 0 ... +5 V ( $R_i = 34\text{ k}\Omega$ )
Analog input 2 (AI2)	resolution 11-bit 0 ... +5 V ( $R_i = 34\text{ k}\Omega$ )
Digital input 1 (DI1)	+3 ... +36 VDC ( $R_i = 12\text{ k}\Omega$ )
Digital input 2 (DI2)	+3 ... +36 VDC ( $R_i = 12\text{ k}\Omega$ )
Digital input 3 (DI3)	+3 ... +36 VDC ( $R_i = 12\text{ k}\Omega$ )

## 4.1.4 Outputs

Digital output 1 (DO1)	open collector, max. 36 VDC (IL < 50 mA)
Digital output 2 (DO2)	open collector, max. 36 VDC (IL < 50 mA)
Digital output 3 (DO3)	push pull, max. 36 VDC (IL < 10 mA)

## 4.1.5 Interfaces

<b>CAN</b>	CAN_H (high); CAN_L (low)	max. 1 Mbit / s (Standard 1 Mbit / s)
<b>RS232</b>	RxD; TxD	max. 115200 bit / s (115200 bit / s)


## 4.1.6 Ambient Conditions

<b>OPERATING TEMPERATURE</b>	-10 ... +45°C
<b>STORAGE TEMPERATURE</b>	-40 ... +85°C
<b>AIR HUMIDITY</b>	Non-condensing 20 ... 80 %
<b>SOUND POWER LEVEL</b>	< 70 dB(A)

## 4.1.7 Connections

<b>POWER SUPPLY</b>	on device: two-row pin connector (2-pole type) Molex Mini Fit Jr plug connector: two-row socket (2-pole type) Molex Mini-Fit Jr. 39-01-2020 crimp contacts: crimp socket Molex Mini-Fit Jr.TM 444-76-1111 (AWG 18-24)
<b>CAN</b>	on device: 3-pole plug JST XARR-03V plug connector: 3-pole socket JST XAP-03V-1
<b>RS232</b>	on device: 4-pole plug JST XARR-04V plug connector: 4-pole socket JST XAP-04V-1
<b>SIGNAL</b>	on device: 12-pole plug JST XARR-12V plug connector: 12-pole socket JST XAP-12V-1

## 4.2 Dimensional Drawing of Pumps






## 4.3 Standard Valve

<b>HOUSING MATERIAL</b>	PEEK
<b>SEALING MATERIAL</b>	FFKM (perfluoroelastomer)
<b>MEDIA TEMPERATURE</b>	0 to +50°C
<b>MAX. VISCOSITY</b>	21 mm <sup>2</sup> /s
<b>INTERNAL VOLUME</b>	approx. 45 µl
<b>MAX. PRESSURE</b>	3 bar
<b>RATED WIDTH</b>	0.6 mm
<b>FLUIDIC CONNECTIONS</b>	¼" – 28 UNF

## 4.4 Dimensional Drawing of Valve



**ATTENTION.** There is a danger of damaging the housing or sealing material. Before using the valve for the first time, please check the chemical compatibility of the media you want to use with the PEEK housing material and the FFKM (perfluoroelastomer) sealing material.

## 4.5 Dosing Performance

The following table provides an overview of the minimum and maximum dosing speeds of the various gear configurations as well as the resulting flow rates, using the example of a 1ml syringe with a 60 mm stroke.

Speeds and flow rates lower than those referred to as pulsation-free, will cause dosing precision to decrease slowly.

GEAR		W/O	14:1	29:1
Min. speed [ $\mu\text{m/s}$ ]		0.065	0.065	0.065
Min. pulsation-free speed [ $\mu\text{m/s}$ ]		14.648	1.042	0.502
Max. speed [ $\text{mm/s}$ ]		89	6.33	3.06
1 ml syringe with 60 mm stroke	Min. flow [ $\mu\text{l/min}$ ]	0.065	0.065	0.065
	Min. pulsation-free flow [ $\mu\text{l/min}$ ]	14.648	1.042	0.502
	Max. flow [ $\text{ml/min}$ ]	89	6.33	3.05

# 5 Transport and Storage

## 5.1 Transport

Please do not lift or transport the modules in assembled state. Transport in assembled condition is only permitted in the original packaging.

Only use the original packaging for transport and shipping.



**ATTENTION.** There is a danger of damaging the device! Do not transport the modules in assembled condition.

## 5.2 Maintenance and Care

When used properly the device is maintenance-free. Should you encounter problems that you cannot fix yourself or which require opening the device, please contact CETONI GmbH to discuss further actions. The device may only be opened by CETONI GmbH or authorized service staff. Violating this rule will void the warranty.

Wipe the device with a moist (not wet) cloth, so that no liquid can drip into the device. In case of heavy soiling you may use a small amount of detergent or alcohol.


# 6 Hardware Operation

## 6.1 Installing a Syringe

The syringe holder on neMESYS modules allows the use of syringes with outside diameters of 6 to 30 mm and a maximum piston stroke of 65 mm.




**IMPORTANT.** Use high-quality glass syringes with outside diameters of 6 to 30 mm to guarantee precise flow rates.



Syringe


Follow this procedure to mount a syringe on the module:

Loosen the knurled screw at the syringe holder. You can now slide the bracket out of the pins and lift it up.



Distance Syringe Holder

Loosen the knurled screw of the piston holder and remove the adapter plate. You can now also lift up the clamping piece.



Distance Piston Holder

To utilize the entire syringe volume, move the piston holder to its maximum forward position through the software. Place an empty syringe onto the remaining base of the syringe holder in such way that the piston touches the piston holder.

The syringe position can be varied somewhat by shifting the piston holder. To do that, loosen the screw with a 4mm Allen wrench.



**ATTENTION.** To avoid damage make sure that the remaining syringe stroke is always the same as or larger than that of the module.



Fixing the Syringe in the Piston Holder


Now reinsert the bracket of the syringe holder, place it on top of the pins and fix the syringe in position using the large knurled screw.



Attachment the Clamping Piece

Put the clamping piece back on and insert the adapter plate matching the piston diameter in such way that the piston head sits between the piston holder and the adapter plate.

The piston is fixed in position by slightly tightening the knurled screw in the back. Make sure that syringe and syringe piston are straight and lined up.



Attachment of the syringe into the Piston Holder




**IMPORTANT.** Syringes, and seals in particular, are wear and tear parts. Please check them on a regular basis and replace them if necessary.

## 6.2 Fluidic / Valve

The neMESYS module can be fitted with an optional valve. This valve allows you to switch the syringe connector between your application (outlet) and a reservoir (intake), thereby enabling automatic refilling of the syringe. You can use the software to adjust the valve in such way that it switches to the inlet automatically during filling.


The three connection points are fitted with a ¼"-28 UNF thread and allow the use of regular HPLC fittings.



Connections

A rocker inside the valve (green) connects the syringe connector (1) to the outlet (2) or the intake (3). In the following image the syringe is connected to the outlet (2), while the intake (red) is closed.

A FFKM membrane (shaded green area) connected to the rocker seals off the fluidic flow. The membrane also limits the valve's operating pressure to 3 bar.



**ATTENTION.** Please observe the maximum operating pressure of 3 bar to avoid damage to the standard valve.



**ATTENTION.** Before using the valve please check the chemical resistance of the PEEK housing material and the FFKM (perfluoroelastomer) sealing material with respect to the fluid you are planning to dose in your application.

## 6.3 Removing Valves


The valve can be removed from the device in a few easy steps. This simplifies the connection of hoses and cleaning, while also enabling you to operate the device without a valve, when not needed.

To remove the valve, push the rocker at the plug and pull it off. When replacing it later on, please make sure to face the plug in the correct direction. The white area should face the valve, while the rocker should point away from the valve.






**ATTENTION**. Please observe the correct placement of the valve plug. The white area should face the valve, while the rocker should point away from the valve.



Rocker


You can then simply lift the valve to remove it.



Disassembly Valve

# 7 Electrical Interfaces

## 7.1 Overview



Interfaces for Connection

All interfaces for connecting the device to a power supply, a PC or an external controller are located on the bottom of the device (figure 1). The following interfaces are available:





Figure 1 – Electrical Interfaces

CONNECTOR	INTERFACES
J1 / J2	Spannungsversorgung 24VDC In / Out
J3 / J4	CAN Schnittstellen In / Out
J5	RS232 Schnittstelle
J6	I/O Signalschnittstelle
J7	Schirm Anschluss

# 7.2 Power Supply (J1 / J2)



## 7.2.1 Pin Assignment


PIN NO.	SIGNAL	DESCRIPTION
1	Ground	Grounding of supply voltage
2	+Vcc	Supply voltage +24 VDC

## 7.2.2 Technical Data

Any power supply can be used, as long as it complies with the following requirements:

POWER SUPPLY REQUIREMENTS	
Output voltage	24 VDC
Output current	Depends on power load and number of modules Typical current draw of one module: 0.3 A Surge current draw of one module: 0.6 A
INFORMATION REGARDING CONTACTS	
Plug connectors	Molex mini-Fit Jr.TM 2-pole type (39-01-2020)
Crimp contacts	Molex mini-Fit Jr.TM crimp socket (444-76-1111)
Crimp tool	Molex crimp tool (69008-0724)

# 7.3 CAN Interface (J3 / J4)



## 7.3.1 Pin Assignment

PIN NO.	SIGNAL	DESCRIPTION
1	CAN_H	CAN high bus line
2	CAN_L	CAN low bus line
3	CAN_GND	CAN Ground

## 7.3.2 Technical Data

Standard type	CAN high-speed, ISO 11898 compatible
Maximum bit rate	1 Mbit/s
Protocol	CANopen DS-301, DS-402
Node ID	Software


## 7.3.3 Connecting neMESYS to CAN bus line CiA DS-102

Pin 1 "CAN_H"	Pin 7 "CAN_H" high bus line
Pin 2 "CAN_L"	Pin 2 "CAN_L" low bus line
Pin 3 "CAN_GND"	Pin 3 "CAN_GND" Masse
Housing "CAN_Shield"	Pin 5 "CAN_Shield" cable shield



**IMPORTANT.** Please observe the maximum Baud rate of your CAN master. The standard Baud rate at the time of shipping is 1 Mbit/s. The CAN bus must be terminated with 2 termination resistors (see section 7.6 on CAN bus termination)

## 7.4 RS-232 Connection (J5)



### 7.4.1 Pin Assignment

PIN NO.	SIGNAL	DESCRIPTION
1	neMESYS TxD	neMESYS RS232 send
2	neMESYS RxD	neMESYS RS232 receive
3	GND	RS232 ground
4	Shield	RS232 shield

### 7.4.2 Technical Data

Maximum input voltage	$\pm 30\text{ V}$
Output voltage	typically $\pm 9\text{ V}$ @ 3k grounded
Maximum bit rate	115200 bit/s (standard 38 400 bit/s)
Internal RS232 driver/receiver	EIA RS232 Standard

### 7.4.3 Connecting neMESYS to a PC

**NEMESYS RS-232 (JST)**

**PC RS-232 (DSUB)**

Pin 1 “neMESYS TxD”

Pin 2 “PC RxD”

Pin 2 “neMESYS RxD”

Pin 3 “PC TxD”


Pin 3 “GND”

Pin 5 “GND”



**IMPORTANT.** Please observe the maximum Baud rate of the RS232 interface on your PC or micro controller. The standard Baud rate setting is 115200 bit/s.

## 7.5 Signal Connection (J6)



“Multi-purpose” digital I/Os and analog inputs are available.

### 7.5.1 Pin Assignment


PIN NR.	SIGNAL	LOG. SIGNAL CONTROL	BESCHREIBUNG
1	AnIN 1	Analog In 1	Analog input 1 (0 – 5V)
2	AnIN 2	Analog In 2	Analog input 2 (0 – 5V)
3	DigIN 1	General Purpose In A	Digital input 1 “multi-purpose”
4	DigIN 2	General Purpose In B	Digital input 2 “multi-purpose”
5	DigIN 3	General Purpose In C	Digital input 3 “multi-purpose”
6	DigOUT 1	General Purpose Out C	Digital output 1 “multi-purpose” / valve voltage
7	DigOUT 2	General Purpose Out D	Digital output 2 “multi-purpose” / valve switching

8	DigOUT 3	General Purpose Out A	Digital output 3 “multi-purpose”
9	GND		Digital signal ground
10	+24 VDC		Auxiliary voltage 24 VDC
11	+5 VDC		Auxiliary voltage output 5 VDC
12	AGND		Analog signal ground

## 7.5.2 Analog Input 1 and 2 (Pin 1 and 2)

As a standard, analog inputs are defined as “general purpose” and can be configured through the software.


AnIN1	Pin [1]
AnIN2	Pin [2]
A_Gnd	Pin [12]
Input voltage range	0 ... 5 VDC
Max. input voltage	0 ... 10 VDC
Input resistance	typisch 34 k $\Omega$ gegen A_GND [12]
A/D converter	11-bit
Resolution	2.49 mV
Bandwidth	250 Hz



## 7.5.3 Digital Inputs 1 and 2 (Pin 3 and 4)

“Multi-purpose” inputs can be freely configured and utilized by the user.

DigIN1	Pin [3] , General Purpose In A
DigIN2	Pin [4] , General Purpose In B
GND	Pin [9]
Type of input	Single ended
Input voltage	0 ... 36VDC
Max. input voltage	-36 ... +36 VDC
Level 0	typically < 0.8 VDC
Level 1	typically > 2.0 VDC
Input resistance	typically 8 k $\Omega$
Input current at level 1	typically 270 $\mu$ A @ 5 VDC
Switching delay	< 300 $\mu$ s



Digital Input 1 and 2


## 7.5.4 Digital Input 3 (Pin 5)

“Multi-purpose” inputs can be freely configured and utilized by the user.

DigIN3	Pin [5], General Purpose In C
--------	-------------------------------



GND	Pin [1]
Eingangsspannung	0 ... 5VDC
Max. Eingangsspannung	-7.5 ... + 12.5 VDC
Pegel 0	typisch < 0.8 VDC
Pegel 1	typisch > 2.0 VDC
Eingangswiderstand	typisch 20 k $\Omega$ gegen GND
Max. Eingangsfrequenz	2.5 MHz




## 7.5.5 Digital Output 1 and 2 (Pin 6 and 7)

### 7.5.5.1 SPEZIFIKATION

Outputs 1 and 2 are used for [switching the valve](#). When the valve is not used the outputs may be free configured and utilized by the user.

DigOUT1	Pin [6], General Purpose Out C
DigOUT2	Pin [7], General Purpose Out D
AGND	Pin [12]
Circuit	Open Collector (internal pull-up resistor 2k2 and diode to +5 VDC)
Switching delay	<3 $\mu$ s




Digital Output 1 and 2

## 7.5.5.2 WIRING EXAMPLES

### DIGOUT „DRAIN“


Maximum input voltage	+36 VDC
Maximum load current	50 mA
Maximum voltage drop	<1.0 V @ 50 mA



Digital Output 1 and 2 Wiring Example „drain“

### DIGOUT „SOURCE“

Output voltage	$U_{out} \approx 5V - 0.75V - (I_{load} \times 2200 \Omega)$
Maximum load current	$I_{load} \leq 2 \text{ mA}$




Digital Output 1 wiring example „source“

## 7.5.6 Digital Output 3

This output is defined as “high speed” and may be freely configured by the user.

### 7.5.6.1 SPECIFICATION

DigOUT3	Pin [8], General Purpose Out A
GND	Pin [9]
Circuit	Push-pull stage
Switching delay	<10 ns



Digital Output 5


### 7.5.6.2 WIRING EXAMPLES

#### DIGOUT „DRAIN“

Maximum input voltage	$U_{in} \leq 3.3 \text{ VDC}$
Maximum load current	$I_{max} \leq 24 \text{ mA}$

Maximum voltage drop

$$U_{\text{drop}} < 0.55 \text{ V} + (I_{\text{max}} \times 100 \Omega)$$




#### DIGOUT „SOURCE“

Output voltage

$$U_{\text{out}} \approx 3.3 \text{ V} - 0.75 \text{ V} - (I_{\text{load}} \times 100 \Omega)$$

Maximal load current

$$I_{\text{load}} \leq 24 \text{ mA}$$







## 7.6 CAN-Bus Termination

The CAN-Bus must be terminated on both ends with a 120 Ω termination resistor.

## 7.7 Switching the Integrated Valve

The integrated valve is switched by [digital outputs 1 and 2](#). Switching of the valve is done through digital output 2. Digital output 1 can be used to lower the voltage during the activated state in order to

prevent the valve from heating up due to the coil. The following logic table shows the relevant valve conditions for all possible signal combinations of digital outputs 1 and 2:

DIG. OUT 1 VALVE VOLTAGE	DIG. OUT 2 VALVE SWITCHING	COIL	LED	VALVE STATE
0	0	aus	aus	
0	1	Spannung abgesenkt	weniger hell	
1	0	aus	aus	
1	1	Spannung 24V	hell	

This means that for switching the valve you set both outputs to 1. This switches the valve with the maximum switching voltage of 24 V. If you want to prevent the valve from heating up, switch digital output 1 to 0 after approximately 1 second. This causes the voltage of the coil current to decrease (15 V) and thereby lessens the heating effect.

To switch the valve to the other direction simply set both outputs to 0.


After turning the device on, both outputs are set to 0.



**IMPORTANT.** Activating the coil in the valve for an extended period of time causes the coil to heat up, which, in case of low flow rates, may warm up the medium flowing through the valve. Lowering the valve voltage after switching lessens this heating-up effect.

# 8 Cable Set

## 8.1 Power Cable




PINS A	SIGNAL	PINS B
-	Ground – Power GND	1
+	Supply voltage +24 VDC	2
C	Shield	C

### TECHNICAL DATA

Cable size	2 x 1.5 mm²
Side A	Cable end sleeves 1.5 mm²
Side B	Molex Mini-Fit Jr. 39-01-2020 Molex Mini-Fit Jr. crimp contact socket 444-76-xxxx

## 8.2 Power Connection Cable




PINS A	SIGNAL	PINS B
1	Ground – Power GND	1
2	Supply voltage +24 VDC	2

## TECHNICAL DATA

Cable Size	2 x 1.3 mm <sup>2</sup>
Side A	Molex Mini-Fit Jr. 39-01-2020 Molex Mini-Fit Jr. crimp contact socket 444-76-xxxx
Side B	Molex Mini-Fit Jr. 39-01-2020 Molex Mini-Fit Jr. crimp contact socket 444-76-xxxx

## 8.3 CAN Cabel (terminated DSub Socket)



PINS A	SIGNAL	PINS B
Housing	CAN shield	C
3	CAN GND	3
2	CAN low	2
7	CAN high	1


## TECHNICAL DATA

Cable size	3 x 0.125 mm <sup>2</sup> , twisted pair, shielded
Side A	D-Sub socket DIN 41652, 9-pole type, with attachment screws, 120Ω resistor between CAN high (7) and CAN low (2)
Side B	3-pole socket JST XAP-03V-1



**IMPORTANT.** The socket contains a 120Ω resistor placed between CAN high and CAN low, acting as a CAN-Bus termination.

# 8.4 CAN Cabel (DSub Plug)




PINS A	SIGNAL	PINS B
Housing	CAN shield	C
3	CAN GND	3
2	CAN low	2
7	CAN high	1

TECHNICAL DATA	
Cable size	3 x 0.125 mm², twisted pair, shielded
Side A	D-Sub plug DIN 41652, 9-pole type, with attachment screws
Side B	3-pole socket JST XAP-03V-1



# 8.5 CAN Connection Cable




PINS A	SIGNAL	PINS B
1	CAN high	1
2	CAN low	2
3	CAN GND	3
C	CAN Schirm	C

## TECHNICAL DATA

Cable size	3 x 0.125 mm², twisted pair, shielded
Side A	3-pole socket JST XAP-03V-1
Side B	3-pole socket JST XAP-03V-1

# 8.6 CAN Terminator




PINS A	SIGNAL	PINS B
1	CAN high	1
2	CAN low	2
3	CAN GND	3

## TECHNICAL DATA

Termination	120Ω resistor between CAN high (1) and CAN low (2)
Plug connector	3-pole socket JST XAP-03V-1

# 8.7 RS232 Cable (DSub Socket)




PINS A	SIGNAL	PINS B
Housing	Shield	4
5	Ground	3
2	PC RxD / neMSYSY TxD	1
3	PC TxD / neMSYSY RxD	2

## TECHNICAL DATA

Cable size	3 x 0.125 mm <sup>2</sup> , twisted pair, shielded
Side A	D-Sub plug DIN 41652, 9-pole type, with attachment screws
Side B	4-pole socket JST XAP-04V-1

# 8.8 USB-to-CAN Adapter



PINS A	SIGNAL
Housing	CAN shield
3	CAN GND
2	CAN low
7	CAN high

TECHNICAL DATA	
Cable Size	3 x 0.125 mm², twisted pair, shielded
Side A	D-Sub plug DIN 41652, 9-pole type, with attachment screws
Side B	USB type A

# 9 Setup & Cable Connection

## 9.1 Mechanical Setup

Use the 4 holes (sunk) at the edges of the module's base plate (see section 4.2 for the dimensional drawing of the pump) to install the module into the front plate of your device. Use M3 bolts for attachment.

## 9.2 Introduction to Cable Connection

Two different cable connections are required for installing neMESYS pumps:

1. Connection of the first module to the power supply and the controller (PLC, PC)
2. Connecting the power supply and data line from one module to the next.

Follow these steps to connect your devices:

## 9.3 Step 1 – Connecting the Power Cable

Connect the [power cable](#) to the plug connector [1](#) of your neMESYS OEM module. Connect the other end of the cable to your power supply (+24 VDC). The required output current is load-dependent (typical current for one module  $\approx 0.3$  A; surge peak current  $\approx 0.6$  A).

## 9.4 Step 2 – Data Connection

### 9.4.1 Introduction

Data connection to the PC can be established through RS232 or CAN bus. The cable connection for both connection types is described in the following sections. Use the section relevant to your type of connection.

## 9.4.2 Data Connection Through CAN

You may choose from two connection options to connect neMESYS modules through CAN:

- *Option 1* – If you use neMESYS pumps exclusively and there is no existing CAN bus or you would like to connect the pumps directly to a PC using CAN, please connect the [CAN cable \(terminated DSub socket\)](#) to CAN interface [J3](#) of the first neMESYS OEM module. Then connect the DSub socket of the cable directly to the respective plug of the [USB-to-CAN adapter](#).
- *Option 2* – If you want to integrate the neMESYS pumps into an existing CAN bus, please use the [CAN cable \(DSub plug\)](#). Connect the cable to CAN interface [J3](#) of the first neMESYS OEM module. Then connect the DSub plug to your existing CAN bus.



**IMPORTANT.** The CAN cable (terminated DSub socket) already contains a 120Ω bus termination resistor to terminate one side of the CAN bus. When using the CAN cable (DSub plug), you have to establish the termination yourself (see CAN terminator)

## 9.4.3 Data Connection Through RS232

Connect the [RS232 cable \(DSub plug\)](#) to RS232 connection [J5](#) of the first neMESYS OEM module. Connect the 9-pole DSub socket at the other end of the cable to the respective plug connector of your controller (PC, PLC).

## 9.5 Step 3 – Power Supply of Additional Modules

Connect the end of the [power connection cable](#) to power connector [J2](#) of the current OEM-module. Connect the other end to power connector [J1](#) of the next module in your setup. In this way, you can supply power to every additional module.

## 9.6 Step 4 – Data Connection to Additional Modules

Use the [CAN connection cable](#) to connect the data bus to additional neMESYS OEM modules. Connect CAN connector [J4](#) of the first module to CAN connector [J3](#) of the next module. The individual elements of a CAN network should be arranged electrically in line. Using the CAN connection cables you can route the internal CAN bus in a line structure from one module to the next.

## 9.7 Step 5 – CAN Bus Termination

CAN bus lines must be terminated on both ends using a 120 Ohm resistor. You can do this very simply with the neMESYS OEM [CAN termination plug](#).

If a neMESYS OEM module is the first device in your CAN network, insert the CAN terminator into CAN connector [J3](#). If you use the [CAN cable \(terminated DSub socket\)](#), this will not be necessary since the resistor is already integrated into the cable's DSub socket. If a neMESYS OEM module is the last module in your CAN network, insert the CAN terminator into CAN connector [J4](#) of this module.

# 10 Disposal

Please send your old devices back to CETONI GmbH. We will take care of proper disposal.

If necessary, please decontaminate the device before sending it back and attach a completed decontamination declaration with your shipment.